

## Targeted Single-Cell DNA Sequencing using the Tapestri System

For more information, consult the Tapestri Single Cell User Guide (PN 3354)

**NOTE** • **Vortex all reagents.**


- **Thaw -20 °C reagents on ice.**
- **Avoid sources of static and pipette slowly and carefully when handling emulsions.**

### Prepare Single Cell Suspension

**IMPORTANT** Provide debris-free cell suspension with > 80% viability. Final cell concentration is 2,000 - 4,000 cells/μl.

- Using **Cell Buffer** dilute cells to ~ 3,000 cells/μl in a total volume of at least 100 μl.

### Encapsulate Cells

- Mount the Base Plate onto the Tapestri Instrument.
- Place 0.2 mL Axygen MAXYmum Recovery PCR tube into the middle of the slot at the left side of the Base Plate.
- Place the DNA cartridge onto the base plate. 
- Pipette **100 μl of Lysis Buffer** (●) into reservoir 1.
- Pipette **100 μl of Cell Suspension (3,000 cells/μl)** into reservoir 2.
- Pipette **200 μl of Encapsulation Oil** into reservoir 3.
- Apply DNA Gasket** on top of the cartridge and close instrument lid.
- Run the Cell Encapsulation program by pressing **Step 1: Encapsulation** on the Tapestri instrument touch screen.
- Once program is completed, press **DONE** and remove cartridge and collection tube from base plate.
- Use a **gel loading tip** to carefully **remove up to 100 μl of the oil layer** at the bottom of the sample leaving no more than 5 μl of oil.

### Lysis + Protease Digest

- Place sample tube into thermal cycler and run the Lysis/Digest protocol:


Step	Temperature	Time
1	50 °C	60 min
2	80 °C	10 min
3	4 °C	HOLD

- Once complete, store the lysed and digested sample at 4 °C until required in the next step.

### Barcode Cells


**IMPORTANT** Protect Barcoding Beads from light.

#### Prime the DNA Cartridge

- Label eight 0.2 mL Axygen MAXYmum Recovery PCR tubes with numbers 1 to 8 and load them into the eight bottom slots of the Base Plate. 
- Mount the DNA Cartridge again onto the Base Plate.
- Pipette **200 μl of Electrode Solution** into reservoirs 4 and **500 μl of Electrode Solution** into reservoirs 5.
- Apply DNA Gasket** on top of the cartridge and close instrument lid.
- Run the Priming program by pressing **Step 2: Barcoding** on the Tapestri Instrument touch screen.
- Once program is completed, press **DONE** and proceed to step 9.
- In the meantime prepare 300 μl Barcoding Master Mix as shown in the table on the right top side of this page.
- Briefly vortex the Barcoding Master Mix and centrifuge to collect the contents and store on ice. *Continued on second column of this page.*

	AML	Myeloid
Reagent	Volume (μl)	Volume (μl)
Barcoding Mix	240	235
Barcoding Mix Additive	60	50
Rv Primer Pool Custom	-	15
<b>Total Volume</b>	<b>300</b>	<b>300</b>

### Load the DNA Cartridge

- Retrieve lysed and digested sample from thermal cycler.
- Pipette the entire **cell encapsulation sample (~80 μl)**, including any oil at the bottom, into reservoir 6. 
- Retrieve Barcoding Beads and vortex tube for 10 seconds at high speed.
- Carefully pipette **200 μl of Barcoding Beads** into reservoir 7.
- Pipette **250 μl of Barcoding Master Mix** into reservoir 8.
- Pipette **1.25 mL of Barcoding Oil** into reservoir 9.
- Apply DNA Gasket** on top of the cartridge and close instrument lid.
- Run the Cell Barcoding program by pressing **NEXT** on the Tapestri Instrument touch screen.
- Once the program is completed, press **DONE** and remove the base plate together with the cartridge and the eight collection tubes
- Remove the DNA Cartridge from the Base Plate.
- Use a **gel loading tip** to carefully **remove up to 90 μl of oil** from the bottom layer of each sample tube.
 

**NOTE** Volumes may vary. The final volume of oil must be ~35 μl per tube & the total volume in each tube must not exceed 100 μl.

### UV Treatment

- Place entire Base Plate with tubes (closed) containing emulsions on ice.
- Place a UV light on top of the ice bucket and expose the samples to non-irradiating **UV light for 8 minutes**.

**IMPORTANT** Make sure to use the **Analytik Jena Blak-Ray XX-15L UV light source**.

### Targeted PCR Amplification

- After UV exposure, remove the Base Plate from the ice, transfer the sample tubes to the thermal cycler, and run the Targeted PCR protocol:

		AML	Myeloid	
Step	Temperature	Time	Time	Cycle
1	95 °C	10 min	10 min	20
2	95 °C	30 sec	30 sec	
3	72 °C	10 sec	10 sec	
4	61 °C	4 min	12 min	
5	72 °C	30 sec	30 sec	
6	72 °C	2 min	2 min	
7	4 °C	HOLD	HOLD	

### Break Emulsions

- Add **5 μl of Extraction Agent** (●) to each tube. Vortex and quick spin.
- Incubate at room temperature for **3 minutes**.
- Add **100 μl of nuclease-free water** to each sample.
- Briefly vortex and spin down for 10 seconds in a mini centrifuge to separate the aqueous and oil layers.
- Pipette 100 μl of the aqueous top layer in each tube into a new 0.2 mL PCR tube. **Do not transfer any oil.**

## Clean Up PCR Product

**NOTE** Equilibrate Ampure XP reagent to room temperature.

Prepare 5 mL fresh 80% ethanol using nuclease-free water.

- For each 100 µl sample, add 75 µl of Ampure XP reagent. Vortex for 10 seconds and quick-spin to collect contents.
  - Incubate tubes at room temperature for 5 minutes.
  - Place on magnet and allow 2 minutes for the beads to separate from solution.
  - Without removing the tubes from the magnet, remove the clear liquid from each tube and discard.
  - To each tube carefully add 200 µl of the freshly prepared 80% ethanol, wait 30 seconds, and remove 200 µl of ethanol without disturbing the SPRI beads.
  - Repeat step 5 once, for a total of two wash cycles.
  - Keeping the tubes on the magnet, remove all residual ethanol from each tube without disturbing the beads.
  - Dry SPRI bead pellets in the tubes on the magnet by incubating at room temperature for 2 minutes. Avoid overdrying beads.
  - Remove the tubes from the magnet. Add 20 µl of nuclease-free water into each tube. Thoroughly resuspend the beads by pipetting up and down several times until no bead clumps are visible.
  - Incubate tubes at room temperature for 2 minutes.
  - Place the tubes onto the magnet and wait for 2 minutes or until solutions are clear.
  - Transfer 15 µl of purified PCR product from each tube to a fresh 0.2 mL PCR tube. Avoid transfer of beads.
  - Store the extracted purified PCR product tubes on ice.
  - Quantify 2 µl of purified PCR product from each sample, using the High Sensitivity Qubit Kit (or equivalent assay) according to the manufacturer's instructions.
- NOTE** The DNA quantity in each sample tube may vary between 0.2 ng/µl to 4.0 ng/µl. If yields are outside this range, contact support@missionbio.com for additional support.
- Store purified PCR product samples at -20 °C until proceeding to the next step.

## Library PCR

- Normalize all eight sample tubes to 0.2 ng/µl in 10 µl each using the Tapestri Sample Quantification Tool (PN 40676).
- In a Pre-PCR area prepare the Library PCR Master Mix as follows:

	Reagent	Volume (µl per 1 reaction)	Volume (µl Pre-Mix)
Pre-Mix	Library Mix	25.0	225.0
	Library P5 Primer (●)	5.0	45.0
	Nuclease-free water	5.0	45.0
	Library P7 Index Primers 1-8 (●)	5.0 per tube	5.0 per tube
	Purified PCR product (0.2 ng/µl)	10.0 per tube	10.0 per tube
	<b>Total Volume</b>	<b>50.0</b>	<b>315.0 (Pre-Mix)</b>

- Aliquot 35.0 µl of Pre-Mix to each of the eight 0.2 mL PCR tubes containing 10 µl of 0.2 ng/µl purified PCR product.
- Add 5 µl of one P7 Index Primer (●) to one sample tube separately. Label the tube with the ID number of the primer added.
- Vortex and quick-spin all eight 50 µl PCR samples to collect contents.
- Transfer samples to a thermal cycler and run the Library PCR protocol:

## Library PCR (continued)

Step	Temperature	Time	Cycle
1	95 °C	3 min	
2	98 °C	20 sec	10
3	62 °C	20 sec	
4	72 °C	45 sec	
5	72 °C	2 min	
6	4 °C	HOLD	

## Clean Up Library PCR Product

**NOTE** Equilibrate Ampure XP reagent to room temperature.

Prepare 5 mL fresh 80% ethanol using nuclease-free water.

- For each 50 µl sample, add 31.5 µl of Ampure XP reagent. Vortex for 10 seconds and quick-spin to collect contents.
- Follow Steps 2 - 13 of Section **Clean Up PCR Product** with the following modifications: Elude the DNA in 12 µl of nuclease-free water (instead of 20 µl) and transfer 10 µl of purified product from each tube to a fresh 0.2 mL PCR tube (instead of 15 µl).

## Quantify and Pool Library

- Quantify all eight libraries using an Agilent Bioanalyzer. DNA 1000 chips may be used with 1 µl of undiluted samples or DNA HS chips may be used with 10x diluted samples.
- Use the Tapestri Sample Quantification Tool (PN 40676) to pool 5 nM of each of the eight libraries in a total volume of up to 40 µl.
- Verify the concentration of the pooled library using a Qubit Fluorometer or equivalent.

## Sequence Tapestri Single-Cell DNA Library (MiSeq Setup)

**IMPORTANT** Illumina's Read 1 primer in well 12 of the MiSeq cartridge needs to be replaced with Mission Bio's Custom Seq Primer (●).

- Dilute 1 µl of 10 nM PhiX stock library in 9 µl dilution buffer (10 mM Tris-HCl with 0.05% Tween 20, pH 8.0) to prepare a 1 nM PhiX library.
- Prepare a fresh aliquot of 0.2 M NaOH by combining 200 µl 1 M NaOH and 800 µl nuclease-free water.
- In two separate 1.5 mL low-bind Eppendorf tubes combine 5 µl of the 5 nM sample library with 5 µl of 0.2 M NaOH or 5 µl of the 1 nM PhiX library with 5 µl of 0.2 M NaOH. Vortex, spin down to collect contents.
- Incubate both libraries for 10 minutes at room temperature.
- Add 990 µl of HT1 Buffer to each library to neutralize the samples. The pooled Tapestri Single-Cell DNA library is now concentrated at 25 pM. The PhiX library is now concentrated at 5 pM.
- Prepare the final sample library as follows:

Reagent	Volume (µl)
Denaturated Pooled Tapestri DNA Library (25 pM)	532
Denaturated PhiX Library (5 pM)	140
HT1 Buffer	28
<b>Total</b>	<b>700</b>

- Load 650 µl of the sample library into well 17 on the MiSeq cartridge (150 PE chemistry).
- Use a pipette to carefully transfer 400 µl of Illumina Read 1 primer mix from well 12 to an 1.5 mL Eppendorf tube.
- Add 4 µl to Mission Bio Custom Seq Primer (●) to the mix, vortex and spin down to collect contents.
- Reload the entire mix into well 12 of the Illumina MiSeq cartridge.
- Sequence the library following manufacturer's instructions.

For technical support visit [www.missionbio.com/support](http://www.missionbio.com/support)  
Email: [support@missionbio.com](mailto:support@missionbio.com) Phone: +1 415 854 0058