Abstract #1663

Analytical Methods to Identify Tumor Heterogeneity and Rare Subclones in Single-Cell DNA Sequencing Data from Targeted Panels

Sombeet Sahu¹, Manimozhi Manivannan¹, Shu Wang¹, Dong Kim¹, Adam Sciambi¹, Nianzhen Li¹, Nigel Beard¹ ¹Mission Bio, South San Francisco, CA, USA

Abstract

Background: With the advancements of single-cell sequencing technologies it is now possible to interrogate thousands of cells in a single experiment. Single-cell RNA-Seq has been available for several years but high-throughput single-cell DNA analysis is in its infancy. Therefore, it is essential to develop new capabilities for assessing genetic variation present in rare cells and to better understand the role that these cells play in the evolution of tumor progression. To address these challenges and enable the characterization of genetic diversity in cancer cell populations, we developed a novel approach to identify mutation signatures which define subclones present in a tumor population.

| **Methods:** Here we present subclone identification method using data generated on the Mission Bio Tapestri $\mathbb R$ Platform and |analyzed by Tapestri analytical workflow. The pipeline steps involve obtaining raw reads from the sequencer, removing adapters, aligning and mapping the reads, calling individual cells, and identifying genetic variants within each cell. After filtering for high quality variants, we then filter for data completeness to ensure only high quality data is used in downstream processing. The variant-cell matrix is then subjected to identification of subclones. Top variants that define the signature of each subclone are also identified. To validate our methodology, we used two different targeted sequencing panels on model systems with known truth mutations. Our pipeline shows the distinct clusters correlating with titration and cell-line ratios. Cluster associated signature mutations were also identified. The pipeline can be used for multi sample analysis with time-series data from diagnosis to relapse or from primary site to metastasis to understand clonal diversity. These data demonstrate the utility of the Tapestri platform, the analytical pipeline, and associated data visualization capability. Our approach addresses key issues of identifying rare subpopulations of cells down to 0.1%, and transforms the ability to accurately characterize clonal heterogeneity in tumor samples. This high-throughput method advances research efforts to improve patient stratification and therapy selection for various cancer indications.



Conflicts of interest: S.S., M.M., S.W., D.K., A.S., N.L., N.B. are employees and shareholders of Mission Bio, Inc.

	C1	C2
EZH2:chr7:148504854:A/AGACTT	Het	Hom
TP53:chr17:7577581:A/G	Het	WT
TP53:chr17:7578211:C/T	Het	WT
Total	10557 (77.79%)	3015 (22.21%)
1% K562	6467 (99.11%)	58 (0.89%)
50% K562	4090 (58.04%)	2957 (41.96%)





Tapest	ri Solutio	n and F	⁻ u	tur	e Wo	ork		Lea	r
e Solutions from Panel	Design to Insights	Catalog Par	als for	Homato	logy and Solid 1	Tumor			0
Reconstruct cell-level mutation	Explore clone distribution with	Catalog I a		Tiemato	logy and Solid 1	lumoi		Poster	Da
			TAPEST	RI SINGLE-C	ELL DNA PANELS			#022	S.
			# genes	# amplicons	Target regions coverage (Kbp)	Panel Uniformity		#923	Se
Minimissionbio Anni Asnesi Brener Annes Orone Annese Astronom Management		AML Panel	19	50	~8.5	>90%			
EXCREMENTATION INCOMPANY EXCREMENTATION INCOMPANY EXCREMENTATION EXCREMENTATION EXCREMENTATION EXCREMENTATION EXCREMENTATION EXCREMENTATION EXCREMENTATION EXCREMENTATION EXCREMENTATION		CLL Panel	34	286	~54.9	>90%		#940	Se
	2 protectioner and an an and an an and an	Myeloid Panel	47	330	~66.0	>80%			
Margin Page Control of the second se		Tumor Hotspot Panel	59	244	~40.0	>80%		#956	Se
Tapestri	Tapestri								
Pipeline	Insights							#982	Se
								#1167	Se

Using machine learning to optimize assays for single cell targeted DNA ssion B sequencing

Analytical Methods to Identify Tumor Heterogeneity and Rare Subclones in Single Cell DNA Sequencing Data from Targeted Panels