

Error Correction in Single-Cell DNA Sequencing: Finding Rare Allele for MRD Clone

Dong Kim¹, Manimozhi Manivannan¹, Sombeet Sahu¹, Shu Wang¹, Saurabh Gulati¹, Nianzhen Li¹, Adam Sciambi¹ and Nigel Beard¹
¹Mission Bio, South San Francisco, CA, USA

Conflicts of interest: D.K., M.M., S.S., S.W., S.G., A.S., N.L., A.P., H.V.,N.B. are employees and shareholders of Mission Bio, Inc.

Abstract

Background


The Tapestri single-cell DNA sequencing platform contains a small number of errors from polymerase incorporations, structure induced template switching, PCR mediated recombination or DNA-damage, and errors from sequencing. These errors can range from 0.5% to 2% depending on the sequencer. To improve variant calling and minimal residual disease (MRD) detection we developed a novel consensus sequence-based method for correcting the errors and reducing the false-positive rates.

Methods


Here we present a model to correct base position errors in Tapestri single-cell DNA analytical workflow. The error correction method involves 2 steps. First, we train the model with cell BAM files from multiple panels. A pileup is generated around a mismatch position with k-1/2 bases on both sides. During training, a batch of normalized pileups and true reference bases are feed into NN-based model. For each batch, error is calculated between predicted bases and true reference bases. Gradients of this error are calculated with respect to weights of each layer. The gradients are then used to update the weight of the model in a back propagation step. Once a validation accuracy reaches, the trained model and base transition matrix are used to predict a likelihood of a base on given position. If prediction is the same as reference, the model corrects observed non-reference base to reference base. To filter out the noisy reads before passing the data to variant caller, we suppress the quality scores of reads having very low coverage.

To validate this method, we used two different targeted panels on a Latin square model system (4 cell line mixtures with 98.4%, 1%, 0.5% and 0.1% dilutions) with known truth mutations. We ran the Tapestri analytical workflow with and without error correction. With the error correction pipeline, we were able to significantly reduce our false positive rates while maintaining our sensitivity.

Sensitivity and specificity on titration experiment




Results on PBMC sample with known truth




a and c) Two clinical samples were processed through analytical pipeline. Frequency of the variants were counted and compared to before and after error correction. Overall error correction resulted in a decrease in the number of observed variants
b and d) The true variants were known from bulk sequencing. 4/5 of the variants showed same sensitivity before and after. There is one variant with low sensitivity.


Why Single-Cell ?




Mission Bio Tapestri Workflow Overview




Error correction workflow using DNNs



Preparing input data for CNN




Generating substitution matrix



a). Substitution rates are calculated by counting number of bases for a given reference base for valid loci.

b). We observed significant variation in the substitution rates between runs and hence fixed matrix would not work. Substitution rates starts plateauing after sub sampling 4M reads

Tapestri Solution



Catalog Panels for Hematology and Solid Tumor

TAPESTRI SINGLE-CELL DNA PANELS					
	# genes	# examples	Target regions coverage	Reads (millions)	
AML Panel	15	50	~95.5	~950	
Cell Panel	28	280	~95.5	~950	
Normal Panel	47	100	~95.5	~950	
Tumor Hetero Panel	50	244	~95.5	~950	

45-GENE TUMOR HETERO PANEL					
ABL1	CD34	CD38	CD45	CD56	
CD71	CD79B	CD81	CD82	CD83	
CD84	CD86	CD88	CD90	CD94	
CD95	CD97	CD99	CD101	CD103	
CD104	CD105	CD106	CD107	CD108	
CD109	CD110	CD111	CD112	CD113	
CD114	CD115	CD116	CD117	CD118	
CD119	CD120	CD121	CD122	CD123	
CD124	CD125	CD126	CD127	CD128	
CD129	CD130	CD131	CD132	CD133	
CD134	CD135	CD136	CD137	CD138	
CD139	CD140	CD141	CD142	CD143	
CD144	CD145	CD146	CD147	CD148	
CD149	CD150	CD151	CD152	CD153	
CD154	CD155	CD156	CD157	CD158	
CD159	CD160	CD161	CD162	CD163	
CD164	CD165	CD166	CD167	CD168	
CD169	CD170	CD171	CD172	CD173	
CD174	CD175	CD176	CD177	CD178	
CD179	CD180	CD181	CD182	CD183	
CD184	CD185	CD186	CD187	CD188	
CD189	CD190	CD191	CD192	CD193	
CD194	CD195	CD196	CD197	CD198	
CD199	CD200	CD201	CD202	CD203	
CD204	CD205	CD206	CD207	CD208	
CD209	CD210	CD211	CD212	CD213	
CD214	CD215	CD216	CD217	CD218	
CD219	CD220	CD221	CD222	CD223	
CD224	CD225	CD226	CD227	CD228	
CD229	CD230	CD231	CD232	CD233	
CD234	CD235	CD236	CD237	CD238	
CD239	CD240	CD241	CD242	CD243	
CD244	CD245	CD246	CD247	CD248	
CD249	CD250	CD251	CD252	CD253	
CD254	CD255	CD256	CD257	CD258	
CD259	CD260	CD261	CD262	CD263	
CD264	CD265	CD266	CD267	CD268	
CD269	CD270	CD271	CD272	CD273	
CD274	CD275	CD276	CD277	CD278	
CD279	CD280	CD281	CD282	CD283	
CD284	CD285	CD286	CD287	CD288	
CD289	CD290	CD291	CD292	CD293	
CD294	CD295	CD296	CD297	CD298	
CD299	CD300	CD301	CD302	CD303	
CD304	CD305	CD306	CD307	CD308	
CD309	CD310	CD311	CD312	CD313	
CD314	CD315	CD316	CD317	CD318	
CD319	CD320	CD321	CD322	CD323	
CD324	CD325	CD326	CD327	CD328	
CD329	CD330	CD331	CD332	CD333	
CD334	CD335	CD336	CD337	CD338	
CD339	CD340	CD341	CD342	CD343	
CD344	CD345	CD346	CD347	CD348	
CD349	CD350	CD351	CD352	CD353	
CD354	CD355	CD356	CD357	CD358	
CD359	CD360	CD361	CD362	CD363	
CD364	CD365	CD366	CD367	CD368	
CD369	CD370	CD371	CD372	CD373	
CD374	CD375	CD376	CD377	CD378	
CD379	CD380	CD381	CD382	CD383	
CD384	CD385	CD386	CD387	CD388	
CD389	CD390	CD391	CD392	CD393	
CD394	CD395	CD396	CD397	CD398	
CD399	CD400	CD401	CD402	CD403	
CD404	CD405	CD406	CD407	CD408	
CD409	CD410	CD411	CD412	CD413	
CD414	CD415	CD416	CD417	CD418	
CD419	CD420	CD421	CD422	CD423	
CD424	CD425	CD426	CD427	CD428	
CD429	CD430	CD431	CD432	CD433	
CD434	CD435	CD436	CD437	CD438	
CD439	CD440	CD441	CD442	CD443	
CD444	CD445	CD446	CD447	CD448	
CD449	CD450	CD451	CD452	CD453	
CD454	CD455	CD456	CD457	CD458	
CD459	CD460	CD461	CD462	CD463	
CD464	CD465	CD466	CD467	CD468	
CD469	CD470	CD471	CD472	CD473	
CD474	CD475	CD476	CD477	CD478	
CD479	CD480	CD481	CD482	CD483	
CD484	CD485	CD486	CD487	CD488	
CD489	CD490	CD491	CD492	CD493	
CD494	CD495	CD496	CD497	CD498	
CD499	CD500	CD501	CD502	CD503	
CD504	CD505	CD506	CD507	CD508	
CD509	CD510	CD511	CD512	CD513	
CD514	CD515	CD516	CD517	CD518	
CD519	CD520	CD521	CD522	CD523	
CD524	CD525	CD526	CD527	CD528	
CD529	CD530	CD531	CD532	CD533	
CD534	CD535	CD536	CD537	CD538	
CD539	CD540	CD541	CD542	CD543	
CD544	CD545	CD546	CD547	CD548	
CD549	CD550	CD551	CD552	CD553	
CD554	CD555				

45-GENE MYELOID PANEL					
ABL1	CD34	CD38	CD45	CD56	
CD71	CD79B	CD81	CD82	CD83	
CD84	CD86	CD88	CD90	CD94	
CD95	CD97	CD99	CD101	CD103	
CD104	CD105	CD106	CD107	CD108	
CD109	CD110	CD111	CD112	CD113	
CD114	CD115	CD116	CD117	CD118	
CD119	CD120	CD121	CD122	CD123	
CD124	CD125	CD126	CD127	CD128	
CD129	CD130	CD131	CD132	CD133	
CD134	CD135	CD136	CD137	CD138	
CD139	CD140	CD141	CD142	CD143	
CD144	CD145	CD146	CD147	CD148	
CD149	CD150	CD151	CD152	CD153	
CD154	CD155	CD156	CD157	CD158	
CD159	CD160	CD161	CD162	CD163	
CD164	CD165	CD166	CD167	CD168	
CD169	CD170	CD171	CD172	CD173	
CD174	CD175	CD176	CD177	CD178	
CD179	CD180	CD181	CD182	CD183	
CD184	CD185	CD186	CD187	CD188	
CD189	CD190	CD191	CD192	CD193	
CD194	CD195	CD196	CD197	CD198	
CD199	CD200	CD201	CD202	CD203	
CD204	CD205	CD206	CD207	CD208	
CD209	CD210	CD211	CD212	CD213	
CD214	CD215	CD216	CD217	CD218	
CD219	CD220	CD221	CD222	CD223	
CD224	CD225	CD226	CD227	CD228	
CD229	CD230	CD231	CD232	CD233	
CD234	CD235	CD236	CD237	CD238	
CD239	CD240	CD241	CD242	CD243	
CD244	CD245	CD246	CD247	CD248	
CD249	CD250	CD251	CD252	CD253	
CD254	CD255	CD256	CD257	CD258	
CD259	CD260	CD261	CD262	CD263	
CD264	CD265	CD266	CD267	CD268	
CD269	CD270	CD271	CD272	CD273	
CD274	CD275	CD276	CD277	CD278	
CD279	CD280	CD281	CD282	CD283	
CD284	CD285	CD286	CD287	CD288	
CD289	CD290	CD291	CD292	CD293	
CD294	CD295	CD296	CD297	CD298	
CD299	CD300	CD301	CD302	CD303	
CD304	CD305	CD306	CD307	CD308	
CD309	CD310	CD311	CD312	CD313	
CD314	CD315	CD316	CD317	CD318	
CD319	CD320	CD321	CD322	CD323	
CD324	CD325	CD326	CD327	CD328	
CD329	CD330	CD331	CD332	CD333	
CD334	CD335	CD336	CD337	CD338	
CD339	CD340	CD341	CD342	CD343	
CD344	CD345	CD346	CD347	CD348	
CD349	CD350	CD351	CD352	CD353	
CD354	CD355	CD356	CD357	CD358	
CD359	CD360	CD361	CD362	CD363	
CD364	CD365	CD366	CD367	CD368	
CD369	CD370	CD371	CD372	CD373	
CD374	CD375	CD376	CD377	CD378	
CD379	CD380	CD381	CD382	CD383	
CD384	CD385	CD386	CD387	CD388	
CD389	CD390	CD391	CD392	CD393	
CD394	CD395	CD396	CD397	CD398	
CD399	CD400	CD401	CD402	CD403	
CD404	CD405	CD406	CD407	CD408	
CD409	CD410	CD411	CD412	CD413	
CD414	CD415	CD416	CD417	CD418	
CD419	CD420	CD421	CD422	CD423	
CD424	CD425	CD426	CD427	CD428	
CD429	CD430	CD431	CD432	CD433	
CD434	CD435	CD436	CD437	CD438	
CD439	CD440	CD441	CD442	CD443	
CD444	CD445	CD446	CD447	CD448	
CD449	CD450	CD451	CD452	CD453	
CD454	CD455	CD456	CD457	CD458	
CD459	CD460	CD461	CD462	CD463	
CD464	CD465	CD466	CD467	CD468	
CD469	CD470	CD471	CD472	CD473	
CD474	CD475	CD476	CD477	CD478	
CD479	CD480	CD481	CD482	CD483	
CD484	CD485	CD486	CD487	CD488	
CD489	CD490	CD491	CD492	CD493	
CD494	CD495	CD496	CD497	CD498	
CD499	CD500	CD501	CD502	CD503	
CD504	CD505	CD506	CD507	CD508	
CD509	CD510	CD511	CD512	CD513	
CD514	CD515	CD516	CD517	CD518	
CD519	CD520	CD521	CD522	CD523	
CD524	CD525	CD526	CD527	CD528	
CD529	CD530	CD531	CD532	CD533	
CD534	CD535	CD536	CD537	CD538	
CD539	CD540	CD541	CD542	CD543	
CD544	CD545	CD546	CD547	CD548	
CD549	CD550	CD551	CD552	CD553	
CD554	CD555				

54-GENE CLL PANEL					
ATM	CD20	CD22	CD23	CD25	
CD38	CD44	CD45	CD48	CD49	
CD50	CD52	CD54	CD56	CD57	
CD58	CD59	CD61	CD62	CD63	
CD64	CD65	CD67	CD68	CD69	
CD70	CD71	CD72	CD73	CD74	
CD75	CD76	CD77	CD78	CD79	
CD80	CD81	CD82	CD83	CD84	
CD85	CD86	CD87	CD88	CD89	
CD90	CD91	CD92	CD93	CD94	
CD95	CD96	CD97	CD98	CD99	
CD100	CD101	CD102	CD103	CD104	
CD105	CD106	CD107	CD108	CD109	
CD110	CD111	CD112	CD113	CD114	
CD115	CD116	CD117	CD118	CD119	
CD120	CD121	CD122	CD123	CD124	
CD125	CD126	CD127	CD128	CD129	
CD130	CD131	CD132	CD133	CD134	
CD135	CD136	CD137	CD138	CD139	
CD140	CD141	CD142	CD143	CD144	
CD145	CD146	CD147	CD148	CD149	
CD150	CD151	CD152	CD153	CD154	
CD155	CD156	CD157	CD158	CD159	
CD160	CD161	CD162	CD163	CD164	
CD165	CD166	CD167	CD168	CD169	
CD170	CD171	CD172	CD173	CD174	
CD175	CD176	CD177	CD178	CD179	
CD180	CD181	CD182	CD183	CD184	
CD185	CD186	CD187	CD188	CD189	
CD190	CD191	CD192	CD193	CD194	
CD195	CD196	CD197	CD198	CD199	
CD200	CD201	CD202	CD203	CD204	
CD205	CD206	CD207	CD208	CD209	
CD210	CD211	CD212	CD213	CD214	
CD215	CD216	CD217	CD218	CD219	
CD220	CD221	CD222	CD223	CD224	
CD225	CD226	CD227	CD228	CD229	
CD230	CD231	CD232	CD233	CD234	
CD235	CD236	CD237	CD238	CD239	
CD240	CD241	CD242			