

Enabling single cell analysis of copy number variation in breast cancer



Jacqueline Marin¹, Saurabh Parikh¹, Arnaud Da Cruz Paula², Shirin Issa Bhaloo², Britta Weigelt², Alex Li¹, Khushali Patel¹, Ania Wronski¹, Anup Parikh¹, Jorge S Reis-Filho²

¹ Mission Bio, South San Francisco, CA, USA. All are employees and shareholders.

² Memorial Sloan Kettering Cancer Center, New York, NY



The Tapestri Solution

Introduction

- There is a lack of data on CNV in single cells in solid tumor samples
- We demonstrate single cell DNA sequencing using the Tapestri instrument for solid tumor samples (Figure 1)
- Informatics pipeline that can detect CNV and SNV in breast cancer tissue (Figure 2)
- Ability to uncover clonal architecture, resistance mechanisms and progression of disease at a single cell level.




Figure 1: Tapestri workflow

Methods

- Panel designed in collaboration with Dr. Jorge Reis-Filho to capture SNV and CNV.
- Minimum of 5 amplicons designed per CNV region
- Nuclei Isolation protocol was optimized to eliminate nuclei aggregates
- 3 patient samples run on Tapestri instrument and analyzed via pipeline (Fig 1 & 2)




Figure 3: SNV calling

Clones are initially identified using the SNV signature. This patient had one clone with the spike-in Raji cells also identified.




Figure 4: Chromosomal changes identified at an arm level

Each dot represents an amplicon assayed using single cell DNA-sequencing. At least 5 amplicons must be changed to call a CNV. This patient sample has a gain of chr1q, chr7q, chr8p and a loss at chr 13q.

Conclusion

Ability to call biologically relevant genomic changes in clinical tissue

Results

- We have developed the only end to end solution for simultaneous detection of SNV and CNV in solid tumors at a single cell level.
- Tapestri was able to identify SNV (Figure 3) and CNV (Figure 4) in clinically relevant samples.
- Copy gain was seen in *EGFR* (Figure 5) which is amplified in 20+% of all breast cancers and is linked to metastasis.
- Copy loss was seen in *BRCA2* and *RB1* (Figure 5). Both events are associated with poor clinical prognosis

Gain	Loss
chr1q	chr13q
chr7q	
chr8p	

Single Cell DNA Sequencing Applications for Breast Cancer

- Clonal Mosaicism and Neoplastic Transition
- Clonal evolution in primary tumors
- Metastatic dissemination & Epithelial to Mesenchymal transition (EMT)
- Drug Resistance Mechanisms

Mission Bio Tapestri Solution

- Simultaneously measures SNVs, indels, and arm-level copy number variations (CNVs) including loss of heterozygosity (LOH) at the single-cell level across up to 10,000 cells
- Unambiguously identifies variant zygosity and mutational co-occurrence
- Detects rare cell populations




Figure 2: Bioinformatics Pipeline




Figure 5: Chromosomal changes identified at a gene level

Each dot represents an amplicon assayed using single cell DNA-sequencing. At least 5 amplicons must be changed to call a CNV. This patient had a gain in *EGFR* and a copy loss in *BRCA1* and *RB1* genes.

Enabling single cell analysis of copy number variation in breast cancer



Jacqueline Marin¹, Saurabh Parikh¹, Arnaud Da Cruz Paula², Shirin Issa Bhaloo², Britta Weigelt², Alex Li¹, Khushali Patel¹, Ania Wronski¹, Anup Parikh¹, Jorge S Reis-Filho²

¹ Mission Bio, South San Francisco, CA, USA. All are employees and shareholders.
² Memorial Sloan Kettering Cancer Center, New York, NY

Introduction

- There is a lack of data on CNV in single cells in solid tumor samples
- We demonstrate single cell DNA sequencing using the Tapestri instrument for solid tumor samples (Figure 1)
- Informatics pipeline that can detect CNV and SNV in breast cancer tissue (Figure 2)
- Ability to uncover clonal architecture, resistance mechanisms and progression of disease at a single cell level.




Figure 1: Tapestri workflow

Methods

- Panel designed in collaboration with Dr. Jorge Reis-Filho to capture SNV and CNV.
- Minimum of 5 amplicons designed per CNV region
- Nuclei Isolation protocol was optimized to eliminate nuclei aggregates
- 3 patient samples run on Tapestri instrument and analyzed via pipeline (Fig 1 & 2)

Tapestri Pipeline Output




Figure 2: Bioinformatics Pipeline

5412 cells, 31 features

Raji 6.4%
 Mixed 3.1%
 Patient 88.5%




Figure 3: SNV calling
 Clones are initially identified using the SNV signature. This patient had one clone with the spike-in Raji cells also identified.

Results

- We have developed the only end to end solution for simultaneous detection of SNV and CNV in solid tumors at a single cell level.
- Tapestri was able to identify SNV (Figure 3) and CNV (Figure 4) in clinically relevant samples.
- Copy gain was seen in *EGFR* (Figure 5) which is amplified in 20+% of all breast cancers and is linked to metastasis.
- Copy loss was seen in *BRCA2* and *RB1* (Figure 5). Both events are associated with poor clinical prognosis




Figure 4: Chromosomal changes identified at an arm level

Each dot represents an amplicon assayed using single cell DNA-sequencing. At least 5 amplicons must be changed to call a CNV. This patient sample has a gain of chr1q, chr7q, chr8p and a loss at chr13q.

Conclusion

Ability to call biologically relevant genomic changes in clinical tissue

Gain	Loss
chr1q	chr13q
chr7q	
chr8p	

Single Cell DNA Sequencing
 Applications for Breast Cancer

- Clonal Mosaicism and Neoplastic Transition
- Clonal evolution in primary tumors
- Metastatic dissemination & Epithelial to Mesenchymal transition (EMT)
- Drug Resistance Mechanisms

Mission Bio
 Tapestri Solution




Figure 5: Chromosomal changes identified at a gene level

Each dot represents an amplicon assayed using single cell DNA-sequencing. At least 5 amplicons must be changed to call a CNV. This patient had a gain in *EGFR* and a copy loss in *BRCA1* and *RB1* genes.

AACR #6635