

Precise Measurement of Viral Copy Number at Single-Cell Resolution for Cell and Gene Therapy Development

Chieh-Yuan (Alex) Li¹, Saurabh Parikh¹, Khushali Patel¹, Qawer Ayaz¹, Yilong Yang¹, Daniel Mendoza¹, Mahir Mohiuddin², Hua-Jun He², John Elliott², Brittany Enzmann¹, Vaneer Pho¹, Sally Floyd¹, Yue Wang¹, Matt Cato¹, Shu Wang¹, Benjamin Schroeder¹
1. Mission Bio, 400 E Jamie Ct, Suite 100, South San Francisco, CA 94080; 2. National Institute of Standards and Technology, Gaithersburg, MD;


Introduction

Cell and gene therapies are transformative solutions for a host of inherited and acquired diseases for which existing interventions are ineffective. Many such therapies rely on the introduction of transgenes into host cells using viral or non-viral vectors. The accurate measurement of gene transfer is critical to the development of therapeutic agents and is a key attribute for assessing their safety and efficacy. Yet, conventional methods for measuring gene transfer lacks the resolution and representation to truly reflect sample composition and either report a population average (bulk) or involve laborious and time-consuming clonal outgrowth which can take weeks. Mission Bio has developed an end-to-end solution from panel design to data analysis for single-cell targeted DNA sequencing. Here, using the Tapestri platform we demonstrate that single-cell DNA sequencing identifies transduced versus non-transduced cells with exceptional accuracy and precision, as well as measure the single-cell level vector copy number (VCN) for populations of thousands of cells with single nucleotide resolution while reducing sample processing time from weeks to days.

Methods


The Tapestri[®] single cell DNA platform utilizes droplet microfluidic technology to rapidly encapsulate, process, and profile up to 20,000 individual cells for multi-analyte detection. The platform is enabled by a novel two-step microfluidic workflow and a high multiplex PCR biochemistry scheme. The two-step microfluidics allows for efficient access to DNA for downstream genomic reactions and provides flexibility to adapt for additional applications and multi-omics (with oligo conjugated antibody during cell preparation). The multiplex PCR chemistry is developed and co-optimized with an AI-powered panel design pipeline and enables direct and efficient amplification of targeted genomic regions within barcoded individual cells. The final products are sequenced on an Illumina sequencing instrument (Figure 1).

Figure 1: Tapestri workflow



Characterized, experimental, lentiviral transduced Jurkat clonal cell lines with vector copy number of 0, 1, 2, 3 and 4 copies were used in this study¹. A custom panel of amplicons was designed to enable the vector copy number analysis. Pure clonal cell lines of VCN 0-4 and three (3) mixtures with non-transduced Raji cell spike-in was processed for single-cell sequencing using Mission Bio's Tapestri platform in triplicates. Tapestri data from a control cell line with known VCN, as well as unknown samples, were combined and analyzed with a prototype copy number caller to produce single-cell level VCN calls, as well as a population average VCN and a percent transduction for each sample. A median of 5,116 cells were analyzed for each Tapestri sample.

Figure 2: Lentiviral transduced VCN cell lines and vector assay design¹




Results

Tapestri vector copy number assessment tool is precise and accurate

Table 1: Combinations of controls and test samples' average VCN

Average copy number for various combinations of controls and test samples	Test Sample											
	VCN1 rep1	VCN1 rep2	VCN1 rep3	VCN2 rep1	VCN2 rep2	VCN2 rep3	VCN3 rep1	VCN3 rep2	VCN3 rep3	VCN4 rep1	VCN4 rep2	VCN4 rep3
VCN1 rep1	1.00	1.00	1.00	2.10	2.00	2.00	3.00	3.00	3.25	4.34	4.46	4.33
VCN1 rep2	1.00	1.00	1.00	2.00	2.00	2.00	3.00	3.00	4.12	4.23	4.01	4.01
VCN1 rep3	1.00	1.00	1.00	2.01	2.00	2.00	3.00	3.00	3.05	4.09	4.32	4.10
VCN2 rep1	1.01	1.03	1.03	2.01	2.00	2.00	3.00	3.00	3.00	4.08	4.27	4.19
VCN2 rep2	1.02	1.03	1.03	2.01	2.00	2.00	3.00	3.00	3.05	4.18	4.35	4.18
VCN2 rep3	1.01	1.01	1.02	2.00	2.00	2.00	3.00	3.00	3.00	4.09	4.38	4.34
VCN3 rep1	1.09	1.08	1.10	2.10	2.08	2.07	3.00	3.11	3.13	4.37	4.64	4.46
VCN3 rep2	1.03	1.03	1.04	2.00	2.00	2.00	3.00	3.00	3.00	4.18	4.45	4.27
VCN3 rep3	1.06	1.06	1.06	2.02	2.01	2.00	3.00	3.00	3.00	4.22	4.42	4.29
VCN4 rep1	1.00	1.01	1.01	2.00	2.00	2.00	2.99	3.00	3.00	4.01	4.34	4.11
VCN4 rep2	1.00	1.00	1.00	1.94	1.98	1.96	2.77	2.96	2.83	3.85	4.04	3.98
VCN4 rep3	1.00	1.00	1.01	2.00	2.00	2.00	2.92	3.00	2.96	4.00	4.31	4.03

Figure 3: Control vs. average VCN



Clonal populations of engineered VCN cells (VCN1-4) were single cell sequenced and processed through the Tapestri VCN caller individually (n=3). We demonstrate the single-cell VCN caller can use any cell population with a known number of VCN's to characterize a test sample's VCN with high accuracy and precision (Table 1). Fig 3 is a graphical representation of Table 1, where each dot represents the average VCN of the test sample (x-axis) using various control samples. Using VCN caller model established from VCN 1 to VCN 4 each as reference points for determining average VCN produces coefficient of variation of 2.76%, 1.62%, 2.48% and 4.04%; and percent error of 2.17%, 0.43%, 0.05% and 5.56% respectively across all combinations across runs.

Tapestri VCN results correlates well with orthogonal ddPCR technology





Figure 4: Tapestri vs. ddPCR average VCN correlation

X-Y scatter plot of the average VCN for each sample called using Tapestri VCN caller vs. ddPCR measurement. Pure cell lines (VCN 0-4) or mixtures of cell lines (Mixes 1-3) were analyzed in triplicate by both Tapestri and ddPCR. The sample's average VCN was then calculate based on single-cell level VCN (y-axis) and plotted against orthogonal ddPCR measurements (Bio Rad QX200, x-axis, n=3). The average copy number of the population estimated using Tapestri correlates well with ddPCR measurements ($R^2=0.99$). The average and standard deviation of replicate measurements (n=3) are shown in the bottom right table.

Transduction efficiency can be detected with high sensitivity and specificity


Figure 5: Transduction efficiency sensitivity and specificity

Tapestri also provides vector amplicon specific measurements that can be used to determine transduction efficiency. A mixture of non-transduced Raji and vector containing Jurkat cells was used for the study. The sample transduction efficiency was calculated based on detection of vector amplicon reads and validated using the assignment from SNPs specific to Raji (negative) and Jurkat (positive) cell lines. The Tapestri VCN can detect transduction efficiency with >99.6% and >99.9% median specificity and sensitivity.



Single-cell level Vector Copy Number distribution characterization


Figure 6: Single-cell VCN distribution per sample



Besides the measurement of average VCN per sample, Tapestri VCN caller calculates VCN distribution of a particular sample at single-cell level. The pie charts show the VCN distribution calculated for each of the three mixtures, Mix1, Mix2, and Mix3 each with ~10% non-vector Raji cell spiked-in. Note that "Mixture 1" (VCN 1-4's respective proportions = 39.4%, 22.8%, 21.3%, 5.4%) and "Mixture 3" (VCN 1-4's respective proportions = 65.6%, 0%, 1.7%, 22.1%) VCN distributions are drastically different despite having similar average population VCN value of 1.64 and 1.61. Besides the measurement of average VCN per sample, Tapestri VCN caller calculates VCN distribution of a particular sample at single-cell level. The pie charts show the VCN distribution calculated for each of the three mixtures, Mix1, Mix2, and Mix3 each with ~10% non-vector Raji cell spiked-in. Note that "Mixture 1" (VCN 1-4's respective proportions = 39.4%, 22.8%, 21.3%, 5.4%) and "Mixture 3" (VCN 1-4's respective proportions = 65.6%, 0%, 1.7%, 22.1%) VCN distributions are drastically different despite having similar average population VCN value of 1.64 and 1.61. Furthermore, Tapestri also provides the average VCN per transduced cell (bottom of each pie chart in red) due to the ability to measure transduction % and VCN distribution.

Single-cell genotype, SNPs, and transgene vector characterization

Figure 9: Heat-map of genotype, VCN and vector sequence



Conclusion

Using the Mission Bio Tapestri single-cell DNA sequencing platform, this study shows a consistent and reliable assay for in-depth quantification of cell and gene therapy transduction efficiency, single-cell vector copy number distribution, and population average (or, VCN per transduced cell). The single-cell level VCN distribution, as well as single-cell vector sequence validation provides unprecedented resolution and insight to assess the potential functional efficacy and safety for CGT products. Together, these characteristics and data demonstrate the potential to accelerate and streamline both the development and release testing of cell and gene therapy products.