




Aim

- The small population of cancerous cells that remain following treatment, known as measurable residual disease (MRD), is the major cause of relapse in acute myeloid leukemia (AML).
- Usually, these refractory cells have gained additional resistance mutations or changed their surface immunophenotypes in ways that preclude detection and phasing by current gold standard flow cytometry or bulk next-generation sequencing assays.
- For this reason, a multiomic single-cell MRD (scMRD) assay could offer a more comprehensive indicator of relapse and the potential for faster response.
- Here, we present a new scMRD assay with a 0.01% limit of detection that provides single-cell clonal architecture and immunophenotyping to not only identify residual leukemia cells, but also identify putative DNA or protein targets for salvage therapy.
- By combining high sensitivity with multiomics, this assay provides researchers with comprehensive and clinically actionable insights into AML MRD.



The scMRD workflow leverages:

- (i) MACs bead enrichment protocol to enrich for CD34+ and/or CD117+ cells.
- (ii) a DNA and protein panel specifically designed for AML MRD diagnosis and treatment^[1, 2, 3]
- (iii) the ability to multiplex up to three patient samples in a single run.

(iv) a new, automated analysis pipeline to evaluate single-cell multiomics output. The pipeline uses each patient's known germline SNP genotypes for demultiplexing samples.

A Multiomic, Single-Cell Measurable Residual Disease (scMRD) Assay For Phasing **DNA Mutations and Surface Immunophenotypes**

Charlie Murphy PhD*, Kathryn Thompson BS*, Lubna Nousheen MSc, Indira Krishnan PhD, Ben Geller BS, Aaron Llanso BS, Todd Druley MD, PhD, Daniel Mendoza PhD, Adam Sciambi PhD Mission Bio. 400 E Jamie Ct, Suite 101, South San Francisco, CA 94080 * Contributed equally

Results

Fig 2: (A) Genes covered by the scMRD DNA panel. (B) The protein AOCs covered by the scMRD protein panel.

(A)

| ASXL1 | FLT3 | MYC | SF3B1 |
|---------|--------|--------|--------|
| BCOR | GATA1 | MYH11 | SMC1A |
| BRAF | GATA2 | NF1 | SRSF2 |
| CALR | IDH1 | NPM1 | STAG2 |
| CBFB | IDH2 | NRAS | TET2 |
| CBL | IL6R* | PHF6 | TP53 |
| CHEK2 | IP6K1* | PPMID | TRPC4* |
| CSF1R | JAK2 | PTPN11 | U2AF1 |
| CYP4F3* | KIT | RAD21 | UBA1* |
| DNMT3A | KMT2A | RUNX1 | WTI |
| ETV6 | KRAS | SETBP1 | ZEB2* |
| EZH2 | MEIS2* | SF3A1* | ZRSR2 |
| | | | |

| Name | ID | Name | ID |
|-------|-----------------------|--------|-----------------------------|
| CD2 | 0367 anti-human CD2 | CD33 | 0052 anti-human CD33 |
| CD3 | 0034 anti-human CD3 | CD34 | 0054 anti-human CD34 |
| CD7 | 0066 anti-human CD7 | CD38 | 0389 anti-human CD38 |
| CD10 | 0062 anti-human CD10 | CD45RA | 0063 anti-human CD45RA |
| CD11b | 0161 anti-human CD11b | CD56 | 0047 anti-human CD56 (NCAM) |
| CD13 | 0364 anti-human CD13 | CD123 | 0064 anti-human CD123 |
| CD14 | 0081 anti-human CD14 | HLA-DR | 0159 anti-human HLA-DR |
| CD19 | 0050 anti-human CD19 | CD117 | anti-human CD117 (A3C6E2) |
| CD22 | 0393 anti-human CD22 | | |

^{*} Germline SNPs for sample demultiplexing rev: The European LeukemiaNet (ELN). World Health Organization (WHO) (

Smith, Takahashi) identified mutations that co-occu provide a fitness advantage leading to disease progression, treatment resistance and


Table 1: Performance of the scMRD assay on samples containing CD34 + or CD117+ cell lines spiked in at 0.01 %. The enrichment ratio is the ratio of the spike in % post enrichment divided by the spike-in % before enrichment.

(B)

| Sample type | Number of runs | Immuno- phenotype | Spike-in % | Number of spike-in cells (average) | Variant specificity (average) | False positive variants (average) | Fold Enrichment (average) |
|------------------------|-------------------|----------------------|------------|--|-------------------------------------|---|---------------------------------|
| Cell line (KG1) | 15 | CD34+ | 0.01% | 11.3 | 99.9 | 0.36 | 27.8 |
| Cell line (HMC-1.2) | 13 | CD117+ | 0.01% | 10.4 | 99.9 | 0.39 | 25.2 |

Fig 3: Enrichment protocol shows the expected change in cell type percentages in healthy bone marrow. Myeloid progenitors increase and T cells decrease. Each point is from a different Tapestri run. Sticky cells are those that express most proteins, and so are likely dead cells.


Healthy bone marrow cells



Post-enrichment

Pre-enrichment

decreased levels of CD117+.



Conclusions

- cells).

References

- 2. Döhner, Hartmut, et al. "Diagnosis and management of AML in adults: 2022 recommendations from an international expert panel on behalf of the ELN." Blood, The Journal of the American Society of Hematology 140.12 (2022): 1345-1377.



Fig 4: Example bioinformatics analysis of Tapestri[®] data from 1 clinical sample that contains 4 known somatic variants, where the corresponding somatic clones express high levels CD34+ and

• The scMRD assay resolved the clonal architecture identifying multiple leukemic clones with co-occurring mutations.

• The integration of genotype and immunophenotypic further enhanced MRD detection by identifying genotype-specific protein expression patterns.

• The assay demonstrates a limit of detection of 0.01%, specificity of $\geq 99\%$ and a false positive rate of ≤ 3 variants per sample.

 Clinical sample genotype and phenotype align with known truth • Samples run on the new V3 chemistry yielded high cell capture rate (26,252

• By combining high sensitivity & specificity with multiomics, the scMRD AML assay offers a potential scalable solution for comprehensive MRD detection that guides therapeutic decision-making.

. Arber, Daniel A., et al. "International Consensus Classification of Myeloid Neoplasms and Acute Leukemias: integrating morphologic, clinical, and genomic data." Blood, The Journal of the American Society of Hematology 140.11 (2022): 1200-1228.

Heuser, Michael, et al. "2021 Update on MRD in acute myeloid leukemia: a consensus document from the European LeukemiaNet MRD Working Party." Blood, The Journal of the American Society of Hematology 138.26 (2021): 2753-2767.