

Decoding the Mosaicism of Genome Editing with Single-cell Multi-omics Analysis

Saurabh Gulaiti¹, Chieh-Yuan (Alex) Li¹, Benjamin Miltz¹, Qawer Ayaz¹, Joanne Nguyen¹, Yilong Yang, Khushali Patel, Brittany Enzmann¹, Vaneé Pho¹, Matt Cato¹, Kathy Tran¹, Edward Szekeres¹, Shu Wang¹, Benjamin Schroeder¹
 1. Mission Bio, 400 E Jamie Ct, Suite 100, South San Francisco, CA 94080


Introduction

Genome editing has emerged as a revolutionary force within the life sciences, wielding transformative potential in applications such as cell and gene therapy development, disease modeling, and functional genomics. Despite the promise of precision of advanced genome editors, editing outcomes remain largely unpredictable. Different cells subjected to the same editing regimen can yield distinct combinations of edits, varying not only across multiple on-target sites but also between on-target and off-target locations. From the perspective of the fundamental biological unit—a single cell—the zygosity disparity (Mono-allelic, Bi-allelic), heterogeneity in variants (homozygous, heterozygous, compound heterozygous), and their functional impact all contribute to the layer of complexity in the mosaicism of editing outcomes. Current genome editing analyses primarily rely on bulk methods, which, though valuable, provide only an average editing efficiency (at the allelic level) of a population. The nuanced cell-to-cell variation of edits remains elusive within these traditional approaches. Here, we present compelling evidence that the Tapestry[®] Genome Editing (GE) Solution offers a breakthrough in the analysis of knockout (KO) and base editing (BE) experiments. We demonstrate the technology's unique single-cell multi-omics capability to furnish intricate details regarding zygosity and the co-occurrence of on- and off-target edits, thereby affording researchers the granularity needed for precise experimental outcomes.

Methods

The Tapestry[®] single cell DNA platform utilizes droplet microfluidic technology to rapidly encapsulate, process, and profile up to 20,000 individual cells for multi-analyte detection. The platform is enabled by a novel two-step microfluidic workflow and a high multiplex PCR biochemistry scheme. The two-step microfluidics allows for efficient access to DNA for downstream genomic reactions and provides flexibility to adapt for additional applications and multi-omics analysis (with oligo conjugated antibody during cell preparation). The multiplex PCR chemistry is developed and co-optimized with an AI-powered panel design pipeline and enables direct and efficient amplification of targeted genomic regions within barcoded individual cells. The final products are sequenced on an Illumina sequencing instrument (Figure 1).

Figure 1: Tapestry[®] workflow



To evaluate the performance of the pipeline for indel calling (for KO applications), we edited Jurkat cells using CRISPR-Cas9 (Synthego), targeting Programmed Cell Death (PDCD1) and T cell receptor α constant (TRAC). Both edited cell pools and isogenic clones were obtained. The on- and off-target edits were verified through bulk NGS and Inference of CRISPR Edits (ICE) analysis. For base editing applications, where base substitutions are the primary edit types, we assess the pipeline's performance on SNV calls using well characterized cell lines, GM12878 & GM24385, from the National Institute of Standards and Technology (NIST). Assays were designed targeting each genome coordinates of interest.

For single-cell multi-omics analysis, a mixture of PBMC and CRISPR-Cas9 edited (TRAC) samples were stained with a 45-plex oligo conjugated antibodies (Total Seq-D, Heme Oncology panel, BioLegend) followed by Tapestry[®] workflow. The gene editing results were analyzed using the Tapestry[®] Genome Editing Solution Protein + DNA Pipeline.

Results

Tapestry[®] GE Pipeline has high level of sensitivity, specificity & reproducibility

A total of 13 and 10 different indels and SNVs were assessed between 24 Tapestry single-cell runs. The performance metric was calculated by comparing the target's expected editing status (from expected truth) to the target's editing status called by the pipeline. The pipeline performance metrics at a sample level uses aggregate counts of events across all cells and targets. High sensitivity (SNV = 98.71% with %CV of 0.51%; Indel = 99.87% with %CV of 1.7%), specificity (SNV = 95.27% with %CV of 0.51%; Indel = 99.93% with %CV of 0.03%) and accuracy (SNV = 95.34% with %CV of 0.51%; Indel = 99.93% with %CV of 0.03%); as well as low false positive rate (SNV=4.73%, Indel=0.07%) and false negative rates (SNV=1.29%, Indel =0.13%) were observed for all samples with high reproducibility and a limit of detection of 0.1% (0.06%-0.09%).

Figure 2: Tapestry GE pipeline performance metrics for SNV and Indel




Table 1: performance metrics

Tapestry Performance	Sensitivity	Specificity	False positive rate	False negative rate	Accuracy	Limit of detection
SNVs	98.71%	95.27%	4.73%	1.29%	95.34%	0.10%
Indels	99.87%	99.93%	0.07%	0.13%	99.93%	0.10%

Accurate measurement of co-occurrence and zygosity of edits at cellular level




Figure 3: Co-occurrence and zygosity of edits

The accuracy of the co-occurrence of edits was assessed using a 100% double-edited PDCD1 and TRAC with one off-target clonal cell line. The clonal cell line (Figure 3, right) has on-target PDCD1 and TRAC edits that are compound heterozygous, bi-allelic (PDCD1: +2/+8bp; TRAC: -2/+2bp), as well as harbors a PDCD1, heterozygous, mono-allelic off-target edits (+2/0bp). The assumption is that all cells would contain the aforementioned co-occurring edits. Indeed, in Figure 3 (left graph), the Tapestry[®] GE pipeline measured the co-occurrence across replicates (n=2) and reported an average of 98.6% of cells containing both TRAC and PDCD1 on-target edits and the predicted off-target edit. Furthermore, it correctly assessed that the zygosity of on-target edits were bi-allelically edited, and the off-target was mono-allelically edited, with the edit indel frequency close to the expected 50:50 range for a diploid cell (46.01%-53.12%, Figure 3, right side, above each allele represents the measured allele frequency). An additional locus that was not targeted for genome editing was included as a negative control (WT). The pipeline also reports the sequence of each indels (not shown).

Sensitive report of indel length distribution and predicted editing activity





Figure 4: Indel length distribution, indel % and predicted KO score

The Tapestry GE pipeline reports the sample's indel % and indel length distribution with high sensitivity. In Figure 4, the left bar plot illustrates the allelic level indel length profiles of a TRAC edited pool measured by Tapestry. The measurements of indel length profile are highly reproducible with Pearson's $r > 0.99$ (log weighted) among Tapestry replicates (Figure 4, top right table). The INDEL% (% alleles with indels) and KO score (defined as the percentage of INDELS that are either longer than 21bp in length or, shorter than 21bp in length but not a multiple of 3bp that results in a coding frameshift) can also be calculated. In this particular edited sample, between replicates, the average INDEL% and inferred KO score are 96.79% (std 1.84%) and 81.60% (std 2.13%), respectively. The indel length, indel %, and KO score calculated by the Tapestry GE pipeline corroborated those calculated by the orthogonal analysis (Synthego).

Direct measurement of editing genotype and functional KO outcome is enabled by single-cell DNA + Protein analysis

The analysis of edited genome and cell-surface protein expression enabled the confirmation of a protein-level knockout in edited cells. Here, we assessed a mixture of TRAC edited Jurkat cells and PBMC with Tapestry DNA + Protein workflow. The Tapestry GE Protein + DNA pipeline reports the sample's edit co-occurrence and zygosity (Figure 3, left), as well as their quantitative surface protein expressions simultaneously. This enables the direct measurement of gene editing outcomes, immunotyping of the samples, and quantitative functional validation of protein expression. To illustrate this, in Figure 5, a mixture of CRISPR-Cas 9, TRAC edited (edit results in Figure 4) Jurkat cells (immortalized human T-cell leukemia cell line, CD3+) are mixed with PBMCs. Indeed, corroborating with previous studies, Tapestry reports TCR α KO disrupts CD3-TCR complex formation and hence impacts surface expression of CD3¹. In addition, from the same data set, more sophisticated single-cell resolution studies, such as per cell allele editing sequence bias (frameshift / indel length) vs. functional protein outcome can be performed.

Figure 5: Tapestry multi-omic analysis validates protein knockout



Conclusion

We have demonstrated that the Tapestry Genome Editing Solution emerges as a robust and comprehensive tool for addressing challenges in genome engineering in decoding the inherent heterogeneity of editing outcomes, both genotypically and functionally. It empowers researchers to navigate the complex landscape of editing variance across thousands of individual cells, with high sensitivity, specificity, accuracy, and reproducibility. Notably, it can detect rare variants down to a remarkable 0.1% of cells, providing unparalleled insight to more sophisticated genome editing endeavors, and their safety and therapeutic efficacy.